


LC-20AB



- | | | |
|---|--|--------------|
| 1 | Suction Filter Assembly | 228-45708-91 |
| 2 | Suction Filter (no tubing) | 228-45707-91 |
| 3 | Inlet Check Valve | 228-45704-91 |
| 4 | Graphite-Filled PTFE Plunger Seal (LC-20AD/AB) | 228-35146-00 |
| 5 | Sapphire Plunger | 228-35601-93 |
| 6 | Diaphragm (2 pack) | 228-32784-91 |
| 7 | Outlet Check Valve | 228-45705-91 |
| 8 | Line Filter | 228-35871-96 |