


# LC-30AD



- |   |   |              |
|---|---|--------------|
| 1 | Suction Filter Assembly                   | 228-45708-91 |
| 2 | Suction Filter (no tubing)                | 228-45707-91 |
| 3 | Inlet Check Valve                         | 228-52964-95 |
| 4 | Polyethylene Plunger Seal and Backup Ring | 228-52711-92 |
| 5 | Sapphire Plunger with Diaphragm           | 228-52069-94 |
| 6 | Diaphragm (2 pack)                        | 228-55272-41 |
| 7 | Outlet Check Valve                        | 228-53334-96 |
| 8 | Line Filter                               | 228-35871-99 |